



Drager DrugCheck 3000

The disposable DrugCheck 3000 test kit is hygienic to use and not likely to be manipulated. The test itself is performed in three easy steps: Swab the oral fluid sample, shake the test kit and start the test after a short waiting period.

- Oral fluid based – non-invasive
- Simultaneous detection of the five most frequent substance classes – two measurement options to detect THC: fast reading or sensitive reading
- Easy to apply
- Pocket sized test



Contact

 1800 506 089

 info@testingdrugs.com.au

 www.testingdrugs.com.au